

ReBorn Essence

Anterior Cervical Plate System

Low Profile

2mm thickness minimizes foreign body sensation

Increased Stability

efficient blockers prevent screw back-out

Anatomical Fit

pre-contoured plate with wide angulation of screws in various footprints allow for ideal fit

ReBorn Essence


Ultra Thin but Secure

The ReBorn Essence® Anterior Cervical Plate System is designed with simplicity, security, and optimal anatomical fit.

The ReBorn Essence® Anterior Cervical Plate System is designed with simplicity, security, and optimal anatomical fit. With 2 mm thickness, the ReBorn Essence Anterior Cervical Fixation System minimizes the sensation of foreign body. The low profile and pre-contoured plate optimizes the anatomical fit. Besides, the intuitive and efficient blocker design provides backout prevention and enhances mechanical stability. The angle-controlled screws in the system come in various footprints and allow flexibility during surgeries.



**Low
Profile**



**Backout
Prevention**



Implant Catalog



■ Cervical Plate

465-26A	D16mm X L26mm
465-28A	D16mm X L28mm
465-30A	D16mm X L30mm
465-32A	D16mm X L32mm
465-34A	D16mm X L34mm
465-36A	D16mm X L36mm
465-38A	D16mm X L38mm
465-40A	D16mm X L40mm
465-42A	D16mm X L42mm
465-44A	D16mm X L44mm
465-46A	D16mm X L46mm
465-48A	D16mm X L48mm
465-50A	D16mm X L50mm
465-52A	D16mm X L52mm
465-54A	D16mm X L54mm
465-58A	D16mm X L58mm
465-62A	D16mm X L62mm
465-66A	D16mm X L66mm
465-70A	D16mm X L70mm
465-74A	D16mm X L74mm
465-78A	D16mm X L78mm
465-82A	D16mm X L82mm
465-84A	D16mm X L84mm
465-86A	D16mm X L86mm
465-90A	D16mm X L90mm

■ Cervical Screw

442-3012A	D3.5mm X L12mm
422-3014A	D3.5mm X L14mm
422-3016A	D3.5mm X L16mm
422-3018A	D3.5mm X L18mm
422-3512A	D4.0mm X L12mm
422-3514A	D4.0mm X L14mm
422-3516A	D4.0mm X L16mm
422-3518A	D4.0mm X L18mm

Instrument Set



1.8mm Awl For Locking Screw
202-0218



2.5mm Hex Screw Driver
202-0825



Plate Holder
202-1301



Plate Bender
202-1201



2.5mm Tap For Locking Plate
202-0525



Handle
202-0101



Temporary Fixation Pin
202-3406



Temporary Fixation Pin Driver
202-3405



Barrier Driver
202-0829



Drill Guide
202-1101



2.5mm Hex Screw Holder
202-0825E



2.3mm Drill Bit
202-1923

Please refer to the use of indication stated in the manual of the ReBorn Essence® Anterior Cervical Fixation System.

BAUI Biotech Co., Ltd.

6F., Sec.1, Zhongxing Rd.,
Wugu Dist., New Taipei City
24872, Taiwan

Tel:+886-2-8976-9538

Fax: +886-2-8976-9530

www.baui.com.tw

Copyright © 2020 BAUI

All rights reserved.

C8804200630

